

Centrifuge Tubes



Your Trusted Partner in Life Sciences



- Abdos centrifuge tubes are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & conforming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- Unique printing technology makes the surface chemical resistant & scratch proof.
- Larger white frosted portion for easy writing.
- Highly transparent, 5% less haze percentage than competitor.
- Hydrodynamic conical bottom for perfect pelleting of sediments.

- Perfect grip on screw cap provides secure control on cap in wet conditions.
- Suitable for temperatures of -80° to +121° C.
- Sterile tubes has been validated according to ENISO 11137-2; SAL 10⁻⁶.



- Even wall thickness to withstand maximum G force, 15000 RCF.
- Specially designed flat thread ensuring leak proof & vacuum tight closure, preventing override of cap.
- Unique plug ensures perfect leak proof seal.



Code	Product Description	Pack Size	No. Per Case
P10401	ABDOS Centrifuge Tubes Conical Bottom, PP/HDPE. 15 ml	25 Pcs / Bag	500
P10402	ABDOS Centrifuge Tubes Conical Bottom, PP/HDPE. 15 ml, Sterile	25 Pcs / Bag	500
P10403	ABDOS Centrifuge Tubes Conical Bottom, PP/HDPE. 50 ml	25 Pcs / Bag	500
P10404	ABDOS Centrifuge Tubes Conical Bottom, PP/HDPE. 50 ml, Sterile	25 Pcs / Bag	500
P10405	ABDOS Centrifuge Tubes Conical Bottom, PP/Amber. 15 ml, Amber	25 Pcs / Bag	500
P10406	ABDOS Centrifuge Tubes Conical Bottom, PP/Amber. 15 ML, Sterile Amber	25 Pcs / Bag	500
P10407	ABDOS Centrifuge Tubes Conical Bottom, PP/Amber. 50 ml, Amber	25 Pcs / Bag	500
P10408	ABDOS Centrifuge Tubes Conical Bottom, PP/Amber. 50 ml, Sterile Amber	25 Pcs / Bag	500
P10419	ABDOS Centrifuge Tubes Conical Bottom, PS/HDPE. 15 ml, Sterile	25 Pcs / Bag	500
P10423	ABDOS Self Standing Centrifuge Tubes, PP/HDPE. 50 ml	25 Pcs / Bag	500
P10424	ABDOS Self Standing Centrifuge Tubes, PP/HDPE. 50 ml, Sterile	25 Pcs / Bag	500
P10425	ABDOS Individually Wrapped Centrifuge Tubes, PP/HDPE. 15 ml, Sterile	1 Ps / 1 Bag	500
P10426	ABDOS Individually Wrapped Centrifuge Tubes, PP/HDPE. 50 ml, Sterile	1 Ps / 1 Bag	500
P10427	ABDOS Individually Wrapped Centrifuge Tubes, PP/HDPE, 50 ml, Sterile, Self Standing	1 Ps / 1 Bag	500

* Abdos Centrifuge tubes are available in environmentally friendly paper rack. For more information Please log on to www.abdoslabware.com/racked-centrifuge-tube-conical-bottom-paper.html .

* For additional product or technical information, please visit www.abdoslabware.com. For distributor information contact labtech@abdosindia.com.

Abdos Labtech Pvt Ltd

720 - 722, DLF Tower A,
Jasola, New Delhi - 110025, India.

T: +91-11-44800980-986

F: +91-11-44800987

E: labtech@abdosindia.com

Europe Office

T: +44 - 7958295643

E: eu.labtech@abdosindia.com