

Technical Data Sheet

Doc No – Ab/TDS/CT-12

Rev No : 03

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Click Together Conical Tube Racks, ABS

Valid for Catalog No. P20218, P20219

1. Description/Specification

1.1 Description:

Click Together Conical Tube Racks, ABS

1.2 Design:

Box with click together design

1.3 Performance Specification:

1.3.1 Capacity: 15 ml x 30 and 50 ml x 20 Centrifuge tube & 15ml x 18tubes, 50ml x 10tubes

1.3.2 Universal Fit for 15 ml and 50 ml Centrifuge Tube Respectively.

1.4 Material:

1.4.1 Acrylonitrile Butadiene Styrene

1.4.2 Comply with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR 177.1020. Olefins polymers.

1.4.3 Complies with the requirement of Indian Standard IS: 10910 on “Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”

1.4.4 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments.

1.5 Color: Yellow/Blue

1.6 Sterilization: Can be autoclaved

1.7 Quality control:

Testing of various characteristic in accordance with the valid Specification



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.8 other Information:

Reusable, Designed for stable use in water bath.

Pack Size: 4 pieces/case

1.9. Basic Features

1.9.1 Made up of Acrylonitrile Butadiene Styrene

1.9.1 Temperature range: -20⁰C to +80⁰C


1.9.3 Autoclavability: Yes, Autoclavable at 15 psi and 121⁰C for 20 minutes

1.9.4 Chemical resistance: See website.

1.9.5 Alphanumeric index aids for the tube identification.

Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121⁰C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com