

# Low Retention Pipette Tips

---



**Your Trusted Partner in Life Sciences**

[www.abdoslabware.com](http://www.abdoslabware.com)




- Abdos last drop low retention pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & conforming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- The universal tip collar fits all popular pipettes.
- Diamond polished moulds for an extremely hydrophobic & smooth surface to ensure lowest retention in its class.
- Advanced super hydrophobic material, confirming retention level is 10 times lesser with viscous samples than standard pipette tips to save an costly reagents.
- The use of the latest injection moulding process produce high precision, thin walled tips with high clarity for better visibility of sample & highly homogeneous surface free of defect sites, preventing non-specific binding.
- Optimal surface for enhanced dispensing behaviour & minimal sample loss of highly viscous liquids or samples containing detergents.
- Sterile tips has been validated according to ENISO 11137-2; SAL 10<sup>-6</sup>.



\* Fig(i)

Abdos Last Drop pipette tips clearly reduce sample loss. As seen in the graph, it clearly indicates that Abdos Tips retained least amount of residue and saves on costly reagents.

- 200 µl Abdos Standard Tips
- 200 µl Abdos Last Drop Tips
- 200 µl Tips of Competitor (A)
- 200 µl Tips of Competitor (T)



\* Fig(i)

Code	Product Description	No. Per Case
P10040	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, Graduated, 0.2 - 10 µl	10 Racks x 96 pcs (960 pcs)
P10041	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, Extended Graduated, 0.2 - 10 µl	10 Racks x 96 pcs (960 pcs)
P10042	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, Bevelled Graduated, 2 - 200 µl	10 Racks x 96 pcs (960 pcs)
P10043	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, Ultra Point Tip Orifice, 2 - 200 µl	10 Racks x 96 pcs (960 pcs)
P100444	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, Graduated, 2 - 300 µl	10 Racks x 96 pcs (960 pcs)
P10045	ABDOS Last Drop Low Retention Pipette Tips, Sterile Racked, 100 - 1000 µl	10 Racks x 96 pcs (960 pcs)

\* Abdos range of Pipette Tips Racked for more details please visit us at [www.abdoslabware.com/low-retention-pipette-tips.html](http://www.abdoslabware.com/low-retention-pipette-tips.html)



Last Drop Low Retention Pipette Tips, Bulk in Resealable Bags



Last Drop Low Retention Pipette Tips, Racked



Last Drop Low Retention Pipette Tips, Sterile Racked



Last Drop Low Retention Pipette Tips, Refill Pack

\* For additional product or technical information, please visit [www.abdoslabware.com](http://www.abdoslabware.com). For distributor information contact [labtech@abdosindia.com](mailto:labtech@abdosindia.com).

## Abdos Labtech Pvt Ltd

720 - 722, DLF Tower A,  
Jasola, New Delhi - 110025, India.

T: +91-11-44800980-986

F: +91-11-44800987

E: [labtech@abdosindia.com](mailto:labtech@abdosindia.com)

### Europe Office

T: +44 - 7958295643

E: [eu.labtech@abdosindia.com](mailto:eu.labtech@abdosindia.com)