

Micro Centrifuge Tubes



Your Trusted Partner in Life Sciences

www.abdoslabware.com



- Abdos micro centrifuge tubes are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & conforming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- Secure seal and precision molded tubes provide a leak proof snap seal closure that allows repeated openings and closings yet resists popping open from internal pressure during centrifugation, boiling and autoclaving.
- Even wall thickness to withstand maximum a force, 20,000 RCF.
- Manufactured using hydrophobic resins on diamond polished moulds resulting in low sample retention.
- Precise lid sealing for lowest MVTR (moist vapour transfer rate) during long-term storage.
- No absorbance @ 260nm & 280nm (safe for genomics & proteomics works).
- Graduation marks that are chemical resistant. providing the researcher with an accurate and easy to read volume measurement.
- Frosted labelling surfaces on the lid and tube.
- Suitable for temprature of -80° to 121° C.



Code

Product Description

No. Per Case


P10201	ABDOS Micro Centrifuge Tubes, PP, 0.5 ml	1000
P10202	ABDOS Micro Centrifuge Tubes, PP, 1.5 ml	500
P10203	ABDOS Micro Centrifuge Tubes, PP, 2.0 ml	500
P10207	ABDOS Micro Centrifuge Tubes, PP, 0.5 ml, Coloured	1000
P10201A	ABDOS Micro Centrifuge Tubes, PP, 0.5 ml, Amber	1000
P10204	ABDOS Micro Centrifuge Tubes, PP, 1.5 ml, Coloured	500
P10204A	ABDOS Micro Centrifuge Tubes, PP, 1.5 ml, Amber	500
P10205	ABDOS Micro Centrifuge Tubes, PP, 2.0 ml, Coloured	500
P10203A	ABDOS Micro Centrifuge Tubes, PP, 2.0 ml, Amber	500

* P10207, P10204 & P10205 is available in Red, Blue, Yellow and Green Colors.

* Abdos range of Micro Centrifuge Tubes Racks for more details please visit us at <http://abdoslabware.com/centrifuge-ware.html>



Rack For Micro Centrifuge Tube



Rack For Micro Centrifuge Tube



Reversible Rack with Cover

* For additional product or technical information, please visit www.abdoslabware.com. For distributor information contact labtech@abdosindia.com.

Abdos Labtech Pvt Ltd

720 - 722, DLF Tower A,
Jasola, New Delhi - 110025, India.

T: +91-11-44800980-986

F: +91-11-44800987

E: labtech@abdosindia.com

Europe Office

T: +44 - 7958295643

E: eu.labtech@abdosindia.com