

Pipette Tips



Your Trusted Partner in Life Sciences

www.abdoslabware.com



- Abdos pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & conforming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- The use of the latest injection moulding processes produces high-precision, thin walled tips with high clarity for better visibility of samples and ensuring perfect seal while reducing the loading and ejection force.
- The universal tip collar fits all popular pipettes.
- Diamond polished moulds for an extremely hydrophobic & smooth surface to ensure lowest retention in its class.
- Perfect air lock grip on pipette by inner radial rings for accurate pipetting.
- The use of multiple cavity tools guarantees high conformity within one batch. The automated moulding and assembling procedures carried out under validated clean-room conditions to eliminate contamination.
- Each wafer of the rack is printed with the effective vol. range and lot no. with product code. It makes easy to keep your work organized and ensure supplies are traceable.
- 5 stage quality control procedures ensure 100% quality & performance.
- Sterile tips has been validated according to ENISO 11137-2; 2012 for SAL 10⁻⁶.



Code	Product Description	No. Per Case
P10117	ABDOS Pipette Tips Sterile Racked, 0.2 - 10 µl, Graduated, Natural	10 Racks x 96 pcs (960 pcs)
P10147	ABDOS Pipette Tips Sterile Racked, 0.2 - 10 µl, Extended Graduated, Natural	10 Racks x 96 pcs (960 pcs)
P10149	ABDOS Pipette Tips Sterile Racked, 2 - 200 µl, Bevelled Graduated, Yellow	10 Racks x 96 pcs (960 pcs)
P10148	ABDOS Pipette Tips Sterile Racked, 2 - 200 µl, Bevelled Graduated, Natural	10 Racks x 96 pcs (960 pcs)
P10111	ABDOS Pipette Tips Sterile Racked, 2 - 200 µl, Ultra Point Tip Orifice, Yellow	10 Racks x 96 pcs (960 pcs)
P10112	ABDOS Pipette Tips Sterile Racked, 2 - 200 µl, Ultra Point Tip Orifice, Natural	10 Racks x 96 pcs (960 pcs)
P10150	ABDOS Pipette Tips Sterile Racked, 2 - 300 µl, Graduated, Natural	10 Racks x 96 pcs (960 pcs)
P10113	ABDOS Pipette Tips Sterile Racked, 100-1000 µl, Blue	10 Racks x 96 pcs (960 pcs)
P10114	ABDOS Pipette Tips Sterile Racked, 100-1000 µl, Natural	10 Racks x 96 pcs (960 pcs)

* Abdos range of Pipette Tips are offered in various configurations for more details please visit us at www.abdoslabware.com/pipette-tips-bulk-pp.html



Pipette Tips, Bulk in Resealable Bags



Pipette Tips, Racked



Pipette Tips, Sterile Racked



Pipette Tips, Refill Pack

* For additional product or technical information, please visit www.abdoslabware.com. For distributor information contact labtech@abdosindia.com.

Abdos Labtech Pvt Ltd

720 - 722, DLF Tower A,
Jasola, New Delhi - 110025, India.
T: +91-11-44800980-986
F: +91-11-44800987
E: labtech@abdosindia.com

Europe Office

T: +44 - 7958295643
E: eu.labtech@abdosindia.com