

Technical Data Sheet

Doc No – Ab/TDS/LH-87

Rev No : 02

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Serological Pipettes Sterile (PS)

Valid for Catalog No. P10601, P10602, P10603, P10604, P10605, P10606

1. Description/Specification

1.1 Description:

Disposable Serological Pipettes Sterile (PS)

1.2 Design:

Transparent Cylindrical Design

1.3 Performance Specifications:

P10601 (1 ml), P10602 (2 ml), P10603 (5 ml), P10604 (10 ml), P10605 (25 ml), P10606 (50 ml)

Accuracy within $\pm 2\%$ At Full volume.

1.4 Material:

1.4.1 Comply with requirement of Directive EU No 10/2011 and requirement of FDA Regulation 21 CFR 177.1640

1.4.2 Materials are manufactured according to the requirements of GMP Regulation (EC) No 2023/2006

1.4.3 Material does not contain any Hazardous chemical substance & compliance with REACH AVHC & RoHS Free.

1.5 Color:

Transparent (Natural)

1.6 Graduations:

Large Cylindrical with temperature, chemical and solvent resistant ink.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com

Technical Data Sheet

1.7 Sterilization:

Sterile catalog number from this lot have been pre-sterilized and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products Requirements for validation and Routine control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation

1.8 Quality control:

Testing of various characteristic in accordance with the valid Specification

1.9 other Information:

Recommended for single use only

2. Basic Features

2.1 Temperature range: -80⁰C to +90⁰C

2.2 Autoclavability: NO

2.3 Chemical resistance: See website.

Autoclaving Guideline for Plastics

Polystyrene can't be autoclaved.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com