

# Technical Data Sheet

Doc No–Ab/TDS/GLP-168

Rev No: 01

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Assorted Workstation Storage Bin, Acrylic**

**Valid for Catalog No. P70104**

### **1. Description/Specification**

**1.1 Description:** Assorted Workstation Storage Bin, Acrylic.

1.1.1 Assorted Workstation Storage Bin can organize most of their commonly used laboratory disposables in one location.

1.1.2 First section of bin is used to store tubes, vials and pipette tips.

1.1.3 Second section of the bin is used to store and dispense tissue paper, Kim wipes and Para film M (supplied with super-sharp serrated blade)

1.1.4 Third section of the bins is used to store and dispense gloves.

### **1.2 Design:**

Workstation Storage Bin

### **1.3 Performance Specifications:**

<b>Code</b>	<b>No. of Compartment</b>	<b>L x W x H (mm)</b>	<b>No. Per Case</b>
<b>P70104</b>	5	255 x 135 x 368	1

### **1.4 Material:**

Poly (methyl methacrylate) (PMMA) Acrylic

**1.5 Color:** Natural

### **1.6 Quality control:**

Testing of various characteristics in accordance with the valid Specification

### **1.7 Other information:**

1.7.1 Strong and lightweight material.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttarakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslaware.com  
**W :** www.abdoslaware.com

# Technical Data Sheet

1.7.2 Good impact strength, higher than both glass & polystyrene.


1.7.3 Transmits up to 92% of visible light (3 mm thickness)

1.7.4 Better chemical stability.

## 1.8. Basic Features

**1.8.1 Autoclavability:** NO, Non autoclavable

**1.8.2 Chemical resistance:** see website.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com