

Technical Data Sheet

Doc No – Ab/TDS/LH-51

Rev No: 00

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Fixed Premium Pipettes, PP

Valid for Catalog No. E10001, E10002, E10003, E10004, E10005, E10006, E10007, E10008, E10009, E10010

1. Description/Specification

1.1 Description: Fixed Premium Pipettes, PP

1.1.1 Abdos fixed premium pipettes are very accurate and gives desired results.

1.1.2 Abdos fixed premium pipettes are lightweight, ergonomic, with low force design.

1.1.3 Volume settings can be clearly read through digital display.

1.1.4 Abdos fixed premium pipettes are easy to calibrate & maintain with the tool supplied.

1.1.5 Abdos fixed premium pipettes are designed to avoid repetitive strain injuries.

1.1.6 Abdos fixed premium pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.

1.1.7 Abdos fixed premium pipettes are suitable for transfer liquid volumes.

1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655/DIN 12650. The quality control involves gravimetric testing of each pipette with distilled water at 22°C.

1.2 Design:

Light weight & simple to use operation.

1.3 Performance Specifications:

Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Sealing & ejecting tips

1.4 Material: PP (polypropylene)



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com

Technical Data Sheet

1.5 Color: White


1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Chemical resistance: see website

1.7. Pipettes

Code	Capacity	No. Per Case
E10001	5 μ l	1
E10002	10 μ l	1
E10003	20 μ l	1
E10004	25 μ l	1
E10005	50 μ l	1
E10006	100 μ l	1
E10007	200 μ l	1
E10008	250 μ l	1
E10009	500 μ l	1
E10010	1000 μ l	1



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata - 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com