

# Technical Data Sheet

Doc No–Ab/TDS/MF-95

Rev No: 01

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Microscope Cover Glass**

**Valid for Catalog No. G10107, G10108, G10109, G10110, G10111, G10112, G10113, G10114**

### **1. Description/Specification**

#### **1.1 Description:** Microscope Cover Glass

- 1.1.1 Microscope Cover glass is made from super white glass sheets for best performance.
- 1.1.2 Microscope Cover Glass Slides are uniform in flatness and sizes for best performance.
- 1.1.3 Designed for manual workflow in histology, cytology, reanalysis and microbiology.
- 1.1.4 Free of worn- out edges, oily films and figure prints.
- 1.1.5 Available in aluminum tropical pack.

#### **1.2 Design:**

Rectangular and square shape design

#### **1.3 Performance Specifications:**

<b>Code</b>	<b>Type</b>	<b>Thickness (mm)</b>	<b>L x W (mm)</b>	<b>No. Per Pack</b>
<b>G10107</b>	Square	1.00 (0.13 – 0.16)	18 x 18	1000
<b>G10108</b>	Square	1.00 (0.13 – 0.16)	20 x 20	1000
<b>G10109</b>	Square	1.00 (0.13 – 0.16)	22 x 22	1000
<b>G10110</b>	Square	1.00 (0.13 – 0.16)	24 x 24	1000
<b>G10111</b>	Rectangle	1.00 (0.13 – 0.16)	22 x 32	1000
<b>G10112</b>	Rectangle	1.00 (0.13 – 0.16)	22 x 40	1000
<b>G10113</b>	Rectangle	1.00 (0.13 – 0.16)	22 x 50	1000
<b>G10114</b>	Rectangle	1.00 (0.13 – 0.16)	22 x 60	1000



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttarakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslalware.com  
**W :** www.abdoslalware.com

# Technical Data Sheet

## 1.4 Material:

Thin and flat glass material.

**1.5 Color:** Transparent

## 1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification


## 1.7 Other information

**1.7.1** ISO 8037/1 compliant.

**1.7.2** Available in aluminum tropical pack.

**1.7.3 Autoclavability:** NO, Non autoclavable

**1.7.4 Chemical resistance:** see website



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com