

Technical Data Sheet

Doc No–Ab/TDS/MF-94

Rev No: 01

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Microscope Slides, Glass

Valid for Catalog No. G10101, G10106

1. Description/Specification

1.1 Description: Microscope Slides, Glass

1.1.1 Microscope slides are made from very uniform, thin and flat glass sheets resistant to corrosion and fogging.

1.1.2 Designed for routine experiments in histology, pathology, hematology, cytology and microbiology for blood, body liquid tests, frozen tissues section examination, biological specimen preparations.

1.1.3 Frosted Microscope slides are double side frosted and can be marked with conventional labeling system, pencil or mark pens.

1.1.4 Frosted area of approx. 20mm depth, resistant to all common chemicals and routine stains used in laboratory.

1.2 Design: Rectangular design, thin and flat glass sheets.

1.3 Performance Specifications:

Code	Type	Thickness (mm)	L x W (mm)	Edges	No. Per Pack
G10101	Plain , Interleaved	1.00 ± 0.2	25 x 75	Ground Edges, 90° Corner	50
G10106	Double Frosted, Interleaved	1.00 ± 0.2	25 x 75	Ground Edges, 90° Corner	50



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.4 Material:

Thin and flat glass material.

1.5 Color: Transparent

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification


1.7 Other information

1.7.1 ISO 8037/1 compliant.

1.7.2 Available with tissue paper interleaved.

1.7.3 Autoclavability: autoclavable

1.7.4 Chemical resistance: see website



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com