

# Technical Data Sheet

Doc No– Ab/TDS/LH-55

Rev No: 02

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Pipette Rack Stand, Acrylic**

**Valid for Catalog No. P70201**

### **1. Description/Specification**

**1.1 Description:** Pipettes Rack Stand, Acrylic

1.1.1 Keep your pipettes within easy inside this convenient rack.

1.1.2 Holds five standard size pipettes.

1.1.3 Racks are mounted on their own stand.

1.1.4 Keep working area or desk free from clutter

### **1.2 Design:**

Pipette Rack Stand

### **1.3 Performance Specifications:**

<b>Code</b>	<b>Places</b>	<b>L x W x H (mm)</b>	<b>No. per Case</b>
<b>P70201</b>	5	190x135x265	1

### **1.4 Material:**

Acrylic, Poly (methyl methacrylate) PMMA

### **1.5 Color:**

Transparent



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com

# Technical Data Sheet

## 1.6 Quality control:

Testing of various characteristic in accordance with the valid Specification

## 1.7. Basic Features

1.7.1 Made up of Acrylic, Poly (methyl methacrylate) PMMA


1.7.2 Strong & Lightweight material.

1.7.3 Good impact strength, higher than both glass and polystyrene.

## 1.8. Basic Features

1.8.1 Autoclavability: NO, Non autoclavable

1.8.2 Chemical resistance: see website.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata - 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com