

# Technical Data Sheet

Doc No – Ab/TDS/CT-03

Rev No : 03

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Rack for Micro Centrifuge Tube, RPP**

**Valid for Catalog No. P20205**

### **1. Description/Specification**

#### **1.1 Description:**

Reusable Rack for Micro Centrifuge Tube, RPP

#### **1.2 Design:**

Compatible Rack for 1.5 ml and 2.0 ml Micro Centrifuge Tube

#### **1.3 Performance Specifications:**

P20205 (1.5ml & 2ml)

#### **1.4 Material:**

1.4.1 Comply with requirement of Directive EU No 10/2011 and requirement of FDA Regulation 21 CFR 177.1520 (c) 3.3b Polymers.

1.4.2 Materials are manufactured according to the requirements of GMP Regulation (EC) No 2023/2006

1.4.3 Material does not contain any Hazardous chemical substance & compliance with REACH AVHC & RoHS Free.

#### **1.5 Color:** White

#### **1.6 Graduations:**

Suitable for 1.5 & 2.0 ml Micro centrifuge Tube

#### **1.7 Sterilization:**

Can be autoclaved at 15 Psi, 121 °C for 20 minutes

#### **1.8 Quality control:**

Testing of various characteristic in accordance with the valid Specification

#### **1.9 other Information:**

Reusable, can be autoclaved up to 10 autoclave cycle.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttarakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com

# Technical Data Sheet

## 1.9 Basic Features


### 1.9.1 Temperature range:

**1.9.2 Autoclavability:** Yes, Autoclavable at 15 psi and 121<sup>0</sup>C for 20 minutes

**1.9.3 Chemical resistance:** See website.

## Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121<sup>0</sup>C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50<sup>0</sup>C before removing product.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com