

# Technical Data Sheet

Doc No Ab/TDS/Eqpt-202

Rev No : 03

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

### **WAVEX – LCD ROTATING MIXER**

**Valid for Catalog No. E11280**

#### **1. Description/Specification**

##### **1.1 Description:**

Wavex Tube Rotator rotating Mixer gives gentle but effective mixing of biological samples to hold in 1.5 ml to 50 ml tubes used in variety of application including immune precipitations, Prevention of blood coagulation and latex diagnostics.

##### **1.2 Design:**

Rotisserie Type

##### **1.3 Performance Specifications:**

- |                                |                              |
|--------------------------------|------------------------------|
| 1.3.1 Voltage (AC)             | : 220~240 Volt               |
| 1.3.2 Frequency (Hz)           | : 50-60                      |
| 1.3.3 Power (W)                | : 40                         |
| 1.3.4 Angle Range              | : 0- 90°                     |
| 1.3.5 Motor type               | : DC Motor                   |
| 1.3.6 Speed range (RPM)        | : 0 ~ 70                     |
| 1.3.7 Speed display            | : LCD                        |
| 1.3.8 Timer range (Min)        | : 1 – 1199                   |
| 1.3.9 Timer display            | : LCD                        |
| 1.4.0 Run type                 | : Timer/Continuous operation |
| 1.4.1 Overall dimension (mm)   | : 220 x 510 x 260 mm         |
| 1.4.2 Weight                   | : 8.0 Kg                     |
| 1.4.3 Permissible ambient Temp | : 5- 40°C                    |



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com

# Technical Data Sheet

1.4.4 Permissible Relative humidity: 80 %

1.4.5 Protection class according to DIN EN60529 : IP21

**1.4 Material:** Body – ABS, Stainless Steel

**1.5 Color:** Dark Blue

## **1.6. Spare Accessories:**

E11281 Rotisserie accessory for 1.5 ml x 48 micro centrifuge tubes horizontal

E11282 Rotisserie accessory for 50 ml x 24 centrifuge tubes horizontally


E11283 Rotisserie accessory for 15 ml x 24 centrifuge tubes horizontally


E11284 Rotisserie accessory for 1.5 ml x 32 micro centrifuge tubes vertically

E11285 Rotisserie accessory for 50 ml x 16 centrifuge tubes vertically

E11286 Rotisserie accessory for 15 ml x 16 centrifuge tubes vertically

E11274 Disk support rod used with disk accessories

  
Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttarakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslaware.com  
**W :** www.abdoslaware.com