

Technical Data Sheet

Doc No – Ab/TDS/RB-47

Rev No: 01

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Weighing Boats, PS

Valid for Cat. No. P70501, P70502, P70503, P70504, P70505, P70506

1. Description/Specification:

1.1 Description: Weighing Boats, PS

1.1.1 Weighing Boats are flexible are with smooth surfaces to provide accurate pour-outs with minimal loss.

1.1.2 Biologically inert and resistant to dilute and weak acids, aqueous solution, alcohols & bases.

1.1.3 These boats can withstand temperature between 10°C to 70°C.

1.1.4 Can be used with microwave moisture analyzes.

1.2 Design:

1.2.1 Weighing Boats are shallow with a wide, flat bottom to resist tipping. They have rounded corners to simplify transfer.

1.3 Performance Specifications:

1.3.1 Weighing Boats Standard, PS

Code	Size	Colour	L x W x H (mm)	No. per Case
P70501	Small	Blue	46x46x8	250
P70502	Medium	Blue	85x85x24	250
P70505	Large	Blue	140x140x22	250

1.3.2 Weighing Boat with Pour Support, PS

Code	Size	Colour	L x W x H (mm)	No. per Case
P70503	Small	White	43x58x13	250
P70504	Medium	White	83x132x26	250
P70506	Large	White	108x183x26	250



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.4 Material:

1.4.1 Complies with requirement of European Pharmacopoeia, 7th Edition (2011) & its supplements 7.5 (07/2012) monograph 3.2.2 Plastic container and closures for Pharmaceutical use.

1.4.2 The product has Passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.5 Color: Blue & White

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification

1.7 Basic Features

1.7.1 Temperature range: 10⁰C to +70⁰C

1.7.2 Autoclavable : No

1.7.3 Chemical resistance: See website.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com