

Technical Data Sheet

Doc No–Ab/TDS/UTL-208

Rev No: 01 Date 04/2023-24

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Abdos Brand – Derma Gaurd Lite Blue Nitrile Gloves – 9.5 inches length

Valid for Catalog No. U20351, U20352, U20353

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Abdos Brand – Derma Gaurd Lite Blue Nitrile Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect your exposure to bacteria, viruses, and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Its unique blue color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.

1.1.3 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.4 Custom design enhances comfort & fit.

1.1.5 These gloves are powder free, minimizing the potential for powder related complications.

1.1.6 Conforms to ASTM D6319, EN 374, EN 455 Standards, EU standard directive EU 2017/745 Class I, GMP, BIS, IMDR 2017, USFDA 510K, ISO 10993, HACCP

1.1.7 Palm Thickness :0.06mm (2.36mil) & Finger Thickness: 0.08mm (3.15mil)

1.1.8 Average Length: 240mm

1.2 Design: Hand gloves design

1.3 Performance Specifications:

Code	Size	No. Per Case
U20351	Small	100
U20352	Medium	100
U20353	Large	100

1.4 Material: Nitrile gloves

1.5 Color: Blue



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.6 Quality control:

Product are manufactured, inspected, and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 These gloves are used for working on lab, machine, and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves

1.7.5 Suitable for food handling

1.7.6 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.7 Manufactured in Aseptic Environment

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.


WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure.

Consult local fire codes prior to storage of flammable liquids.



Technical Data Sheet

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903

Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7