

One of the World's Lightest Pipettes

Elegant Starter Pack

Abdos introduces a Starter Pack for Elegant Micropipettes. These micropipettes are fully autoclavable and feature an ergonomically designed body, offering high accuracy and precision. They are ultralightweight, requiring minimal pipetting force. The soft volume adjustment and piston movement ensure effortless and fatigue-free pipetting, helping to prevent repetitive strain injuries. Each elegant pipette is supplied with an individual test certificate and calibrated in accordance with ISO8655 standards.

**NEW
OFFER**

Micropipettes

Pipette Stand

Pen and Pencil

Gen2 Tip Boxes

Bench-top
Pipette Hangers



Key Features

- ↗ Light & smooth volume aspiration and dispensing
- ↗ One key volume lock for pull and push
- ↗ 4-digit display clearly reads volume setting and provides better incremental value
- ↗ The filter element at tip cone to ensure contamination-free pipetting
- ↗ Chemical and UV-resistant body that allows usability with a variety of samples under UV cabinets
- ↗ Lightweight, ergonomic, low force design to help avoid repetitive strain injuries

Ordering Information

Product Code	Description	Inclusions
E11159	Elegant Starter Pack 01	E11151 (0.1-10µl Pipette) E11153 (10-100µl Pipette) E11156 (100-1000µl Pipette) P10139 (0.1-10µl Gen2 Tip Box) P10140 (1-200µl Gen2 Tip Box) P10109 (100-1000µl Gen2 Tip Box) P70201 (5 Place Linear Pipette Stand) Pipette Bench-top Hangers Abdos Biodegradable Pen & Pencil Set
E11160	Elegant Starter Pack 02	E11152 (2-20µl Pipette) E11154 (20-200µl Pipette) E11156 (100-1000µl Pipette) P10139 (0.1-10µl Gen2 Tip Box) P10140 (1-200µl Gen2 Tip Box) P10109 (100-1000µl Gen2 Tip Box) P70201 (5 Place Linear Pipette Stand) Pipette Bench-top Hangers Abdos Biodegradable Pen & Pencil Set

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships