


# Erlenmeyer Flasks, PC & PETG

Perfect Fit for  
Suspension Cultures



## Features

- Plain or Baffled bottom to improve aeration or reduce shear stress
- HDPE Cap Option of Sealed and Vent for adequate gas exchange
- 0.2µm PTFE hydrophobic membrane for sterile gas exchange and an effective barrier for contamination
- Moulded graduation for accuracy and quick volume check
- Volume options – 125/250/500/1000ml
- Individually wrapped and sterilised by gamma irradiation SAL10-6 as per ISO 11137 guidelines




# Erlenmeyer Flasks, PC & PETG

## Introduction

The Erlenmeyer Flask is a cell culture container, suitable for suspension culture of cell lines with high oxygen requirement (animal cells, bacteria, fungi, plant callus, etc.)


Abdos Life Sciences Erlenmeyer Flasks are available in Polycarbonate (PC) and Polyethylene Terephthalate Glycol (PETG) material, ideal for all shaker & suspension cell culture applications, media preparation, mixing, liquid handling, and storage



## Abdos Erlenmeyer Flask Range



## Quality & Compliance



# Erlenmeyer Flasks Available in Various Sizes



**Baffled: 125/250/500/1000mL**



**Flat: 125/250/500/1000mL**

## Ordering Information


### Erlenmeyer Flasks, PC

Product Code	Capacity	Cap Style	Base Style	Qty/Case
P21200	125ml	Sealed	Flat	24
P21201	125ml	Vent	Flat	24
P21202	250ml	Sealed	Flat	12
P21203	250ml	Vent	Flat	12
P21204	500ml	Sealed	Flat	12
P21205	500ml	Vent	Flat	12
P21206	1000ml	Sealed	Flat	12
P21207	1000ml	Vent	Flat	12
P21212	125ml	Sealed	Baffled	24
P21213	125ml	Vent	Baffled	24
P21214	250ml	Sealed	Baffled	12
P21215	250ml	Vent	Baffled	12
P21216	500ml	Sealed	Baffled	12
P21217	500ml	Vent	Baffled	12
P21218	1000ml	Sealed	Baffled	12
P21219	1000ml	Vent	Baffled	12

### Erlenmeyer Flasks, PETG

Product Code	Capacity	Cap Style	Base Style	Qty/Case
P21220	125ml	Sealed	Flat	24
P21221	125ml	Vent	Flat	24
P21222	250ml	Sealed	Flat	12
P21223	250ml	Vent	Flat	12
P21224	500ml	Sealed	Flat	12
P21225	500ml	Vent	Flat	12
P21226	1000ml	Sealed	Flat	12
P21227	1000ml	Vent	Flat	12


# Other Items of Interest in Cell Culture Range



Cell Culture Flasks




Tissue Culture Plates




Cell Culture Dishes



Biofill Bottle Top  
Filtration System



Serological  
Pipettes



Cell  
Scraper

## What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

📞 011-41081897

✉️ [labtech@abdosindia.com](mailto:labtech@abdosindia.com)

🌐 [www.abdoslifesciences.com](http://www.abdoslifesciences.com)