

## Non-printed and non-assembled

### Non- Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0223	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	5565
PI0227	1.5 ml Self-standing	10.2 mm	44 mm	1000	5565
PI0231	1.5 ml conical	10.2 mm	43.1 mm	1000	5565
PI0235	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	5565
PI0239	2.0 ml conical	10.2 mm	44.3 mm	1000	5565
PI0223A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	6227
PI0223A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	6169
PI0225	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	6169
PI0260A	1.5 ml amber Self-standing	10.2 mm	44 mm	1000	6227
PI0260A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	6169
PI0229	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	6169
PI0270A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	6227
PI0232A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	6169
PI0233	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	6169

\*\* PI0225, PI0229, PI0233, PI0237, PI0241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0272A	2.0 ml amber self-stand	10.2 mm	45.1 mm	1000	6227
PI0236A1	2.0 ml transparent amber self-stand	10.2 mm	45.1 mm	1000	6169
PI0237	2.0 ml coloured self-stand	10.2 mm	45.1 mm	1000	6169
PI0274A	2.0 ml amber conical	10.2 mm	44.3 mm	1000	6227
PI0240A1	2.0 ml transparent amber conical	10.2 mm	44.3 mm	1000	6169
PI0241	2.0 ml coloured conical	10.2 mm	44.3 mm	1000	6169

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0222	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8138
PI0243	1.5 ml Self-standing	10.2 mm	44 mm	1000	8138
PI0244	1.5 ml conical	10.2 mm	43.1 mm	1000	8138
PI0245	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8138
PI0246	2.0 ml conical	10.2 mm	44.3 mm	1000	8138

### Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0224	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6195
PI0228	1.5 ml Self-standing	10.2 mm	44 mm	1000	6195
PI0232	1.5 ml Conical	10.2 mm	43.1 mm	1000	6195
PI0236	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6195
PI0240	2.0 ml Conical	10.2 mm	44.3 mm	1000	6195
PI0224A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	6930
PI0224A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	6930
PI0226	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	6930
PI0228A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	6930
PI0228A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	6865
PI0230	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	6930
PI0232A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	6930
PI0270A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	6865
PI0234	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	6930

\*\* PI0225, PI0229, PI0233, PI0237, PI0241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0236A	2.0 ml amber	10.2 mm	45.1 mm	1000	6930
PI0272A1	2.0 ml transparent amber	10.2 mm	45.1 mm	1000	6865
PI0238	2.0 ml coloured	10.2 mm	45.1 mm	1000	6930
PI0240A	2.0 ml amber	10.2 mm	44.3 mm	1000	6930
PI0274A1	2.0 ml transparent amber	10.2 mm	44.3 mm	1000	6865
PI0242	2.0 ml coloured	10.2 mm	44.3 mm	1000	6930

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring,PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0221	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8715
PI0247	1.5 ml Self-standing	10.2 mm	44 mm	1000	8715
PI0248	1.5 ml conical	10.2 mm	43.1 mm	1000	8715
PI0249	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8715
PI0250	2.0 ml conical	10.2 mm	44.3 mm	1000	8715

## Printed and Non-assembled

### Non- Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30223	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	5933
P30227	1.5 ml Self-standing	10.2 mm	44 mm	1000	5933
P30231	1.5 ml conical	10.2 mm	43.1 mm	1000	5933
P30235	2.0 ml self-standing	10.2 mm	45.1 mm	1000	5933
P30239	2.0 ml conical	10.2 mm	44.3 mm	1000	5933
P30223A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	6594
P30223A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	6536
P30225	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	6536
P30260A	1.5 ml amber Self-standing	10.2 mm	44 mm	1000	6594
P30260A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	6536
P30229	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	6536
P30270A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	6594
P30232A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	6536
P30233	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	6536

\*\*\*P30225, P30229, P30233, P30237, P30241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30272A	2.0 ml amber self-standing	10.2 mm	45.1 mm	1000	6594
P30236A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	6536
P30237	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	6536
P30274A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	6594
P30240A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	6536
P30241	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	6536

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30222	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8505
P30243	1.5 ml Self-standing	10.2 mm	44 mm	1000	8505
P30244	1.5 ml conical	10.2 mm	43.1 mm	1000	8505
P30245	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8505
P30246	2.0 ml conical	10.2 mm	44.3 mm	1000	8505

### Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30224	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6563
P30228	1.5 ml Self-standing	10.2 mm	44 mm	1000	6563
P30232	1.5 ml Conical	10.2 mm	43.1 mm	1000	6563
P30236	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6563
P30240	2.0 ml Conical	10.2 mm	44.3 mm	1000	6563
P30224A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7298
P30224A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7298
P30226	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7298
P30228A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	7298
P30228A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7232
P30230	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	7298
P30232A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7298
P30270A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7232
P30234	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7298

\*\*P30226, P30230, P30234, P30238, P30242 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30236A	2.0 ml amber self-stand	10.2 mm	45.1 mm	1000	7298
P30272A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	7232
P30238	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	7298
P30240A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	7298
P30274A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	7232
P30242	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	7298

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring,PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P30221	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	9083
P30247	1.5 ml Self-standing	10.2 mm	44 mm	1000	9083
P30248	1.5 ml conical	10.2 mm	43.1 mm	1000	9083
P30249	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	9083
P30250	2.0 ml conical	10.2 mm	44.3 mm	1000	9083

## Non-Printed and Assembled

## Non- Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40223	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	5828
P40227	1.5 ml Self-standing	10.2 mm	44 mm	1000	5828
P40231	1.5 ml conical	10.2 mm	43.1 mm	1000	5828
P40235	2.0 ml self-standing	10.2 mm	45.1 mm	1000	5828
P40239	2.0 ml	10.2 mm	44.3 mm	1000	5828
P40223A	0.5 ml amber self-standing	10.2 mm	43.85 mm	1000	6489
P40223A1	0.5 ml transparent amber self-standing	10.2 mm	43.85 mm	1000	6431
P40225	0.5 ml colour self-standing	10.2 mm	43.85 mm	1000	6431
P40260A	1.5 ml self-standing	10.2 mm	44 mm	1000	6489
P40260A1	1.5 ml transparent amber self-standing	10.2 mm	44 mm	1000	6431
P40229	1.5 ml colour self-standing	10.2 mm	44 mm	1000	6431
P40270A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	6489
P40232A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	6431
P40233	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	6431

\*\*P40225, P40229, P40233, P40237, P40241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40272A	2.0 ml amber Self-standing	10.2 mm	45.1 mm	1000	6489
P40236A1	2.0 ml transparent amber Self-standing	10.2 mm	45.1 mm	1000	6431
P40237	2.0 ml coloured Self-standing	10.2 mm	45.1 mm	1000	6431
P40274A	2.0 ml amber conical	10.2 mm	44.3 mm	1000	6489
P40240A1	2.0 ml transparent amber conical	10.2 mm	44.3 mm	1000	6431
P40241	2.0 ml coloured conical	10.2 mm	44.3 mm	1000	6431

## Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40222	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8400
P40243	1.5 ml Self-standing	10.2 mm	44 mm	1000	8400
P40244	1.5 ml conical	10.2 mm	43.1 mm	1000	8400
P40245	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8400
P40246	2.0 ml conical	10.2 mm	44.3 mm	1000	8400

## Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40224	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6458
P40228	1.5 ml Self-standing	10.2 mm	44 mm	1000	6458
P40232	1.5 ml Conical	10.2 mm	43.1 mm	1000	6458
P40236	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6458
P40240	2.0 ml Conical	10.2 mm	44.3 mm	1000	6458
P40224A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7193
P40224A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7193
P40226	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7193
P40228A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	7193
P40228A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7127
P40230	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	6458
P40232A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7193
P40270A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7127
P40234	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7193

\*\*P40226, P40230, P40234, P40238, P40242 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40236A	2.0 ml amber	10.2 mm	45.1 mm	1000	7137
P40272A1	2.0 ml transparent amber	10.2 mm	45.1 mm	1000	7127
P40238	2.0 ml coloured	10.2 mm	45.1 mm	1000	7193
P40240A	2.0 ml amber	10.2 mm	44.3 mm	1000	7193
P40274A1	2.0 ml transparent amber	10.2 mm	44.3 mm	1000	7127
P40242	2.0 ml coloured	10.2 mm	44.3 mm	1000	7193

## Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P40221	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8978
P40247	1.5 ml Self-standing	10.2 mm	44 mm	1000	8978
P40248	1.5 ml conical	10.2 mm	43.1 mm	1000	8978
P40249	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8978
P40250	2.0 ml conical	10.2 mm	44.3 mm	1000	8978

## Printed and Assembled

## Non- Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50223	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6195
P50227	1.5 ml Self-standing	10.2 mm	44 mm	1000	6195
P50231	1.5 ml conical	10.2 mm	43.1 mm	1000	6195
P50235	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6195
P50239	2.0 ml conical	10.2 mm	44.3 mm	1000	6195
P50223A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	6857
P50223A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	6799
P50225	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	6799
P50260A	1.5 ml amber Self-standing	10.2 mm	44 mm	1000	6857
P50260A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	6799
P50229	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	6799
P50270A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	6857
P50232A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	6799
P50233	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	6799

\*\*P50225, P50229, P50233, P50237, P50241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50272A	2.0 ml amber self-standing	10.2 mm	45.1 mm	1000	6857
P50236A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	6799
P50237	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	6799
P50274A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	7495
P50240A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	6475
P50241	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	6799

## Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50222	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	8768
P50243	1.5 ml Self-standing	10.2 mm	44 mm	1000	8768
P50244	1.5 ml conical	10.2 mm	43.1 mm	1000	8768
P50245	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	8768
P50246	2.0 ml conical	10.2 mm	44.3 mm	1000	8768

## Sterile Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50224	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6825
P50228	1.5 ml Self-standing	10.2 mm	44 mm	1000	6825
P50232	1.5 ml Conical	10.2 mm	43.1 mm	1000	6825
P50236	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6825
P50240	2.0 ml Conical	10.2 mm	44.3 mm	1000	6825
P50224A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7560
P50224A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7560
P50226	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7560
P50228A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	7560
P50228A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7495
P50230	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	7560
P50232A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7560
P50270A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7495
P50234	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7560

\*\*P50226, P50230, P50234, P50238, P50242 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50236A	2.0 ml amber self-stand	10.2 mm	45.1 mm	1000	7560
P50272A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	7495
P50238	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	7560
P50240A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	7560
P50274A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	7495
P50242	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	7560

## Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
P50221	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	9345
P50247	1.5 ml Self-standing	10.2 mm	44 mm	1000	9345
P50248	1.5 ml conical	10.2 mm	43.1 mm	1000	9345
P50249	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	9345
P50250	2.0 ml conical	10.2 mm	44.3 mm	1000	9345

## Non-Printed Screw Cap Tube with Looped Cap

### Non-Sterile Screw cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0252	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6400
PI0254	1.5 ml Self-standing	10.2 mm	44 mm	1000	6400
PI0256	1.5 ml conical	10.2 mm	43.1 mm	1000	6400
PI0258	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	6400
PI0260	2.0 ml conical	10.2 mm	44.3 mm	1000	6400
PI0252A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7168
PI0252A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7096
PI0262	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7096
PI0254A	1.5 ml amber Self-standing	10.2 mm	44 mm	1000	7161
PI0254A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7096
PI0264	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	7096
PI0256A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7161
PI0256A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7096
PI0266	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7096

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0258A	2.0 ml amber self-stand	10.2 mm	45.1 mm	1000	7970
PI0258A1	2.0 ml transparent amber self-stand	10.2 mm	45.1 mm	1000	7096
PI0268	2.0 ml coloured self-stand	10.2 mm	45.1 mm	1000	7096
PI0260A	2.0 ml amber conical	10.2 mm	44.3 mm	1000	6227
PI0260A1	2.0 ml transparent amber conical	10.2 mm	44.3 mm	1000	6169
PI0270	2.0 ml coloured conical	10.2 mm	44.3 mm	1000	7096

### Last Drop Low Retention Screw cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0272	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	9359
PI0273	1.5 ml Self-standing	10.2 mm	44 mm	1000	9359
PI0274	1.5 ml conical	10.2 mm	43.1 mm	1000	9359
PI0275	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	9359
PI0276	2.0 ml conical	10.2 mm	44.3 mm	1000	9359

\*\* PI0225, PI0229, PI0233, PI0237, PI0241 are available with red, blue, green, yellow, orange, purple & white color cap options for color coding

### Sterile Screw Cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0253	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	7124
PI0255	1.5 ml Self-standing	10.2 mm	44 mm	1000	7124
PI0257	1.5 ml Conical	10.2 mm	43.1 mm	1000	7124
PI0259	2.0 ml self-standing	10.2 mm	45.1 mm	1000	7124
PI0261	2.0 ml Conical	10.2 mm	44.3 mm	1000	7124
PI0253A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7970
PI0253A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7970
PI0263	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7970
PI0255A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	7970
PI0255A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7896
PI0265	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	7970
PI0257A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7970
PI0257A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7896
PI0267	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7970

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0259A	2.0 ml amber	10.2 mm	45.1 mm	1000	7970
PI0259A1	2.0 ml transparent amber	10.2 mm	45.1 mm	1000	7896
PI0269	2.0 ml coloured	10.2 mm	45.1 mm	1000	7970
PI0261A	2.0 ml amber	10.2 mm	44.3 mm	1000	7970
PI0261A1	2.0 ml transparent amber	10.2 mm	44.3 mm	1000	7896
PI0271	2.0 ml coloured	10.2 mm	44.3 mm	1000	7970

### Last Drop Low Retention Screw Cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0277	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	10022
PI0278	1.5 ml Self-standing	10.2 mm	44 mm	1000	10022
PI0279	1.5 ml conical	10.2 mm	43.1 mm	1000	10022
PI0280	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	10022
PI0281	2.0 ml conical	10.2 mm	44.3 mm	1000	10022

PI0263, PI0265, PI0267, PI0269 & PI0271 are available with red, blue, green, yellow, orange, purple & white color looped cap options for color coding

## Printed Screw Cap Tube with Looped Cap

### Non-Sterile Screw cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0340	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	6823
PI0342	1.5 ml Self-standing	10.2 mm	44 mm	1000	6823
PI0344	1.5 ml conical	10.2 mm	43.1 mm	1000	6823
PI0346	2.0 ml self-standing	10.2 mm	45.1 mm	1000	6823
PI0348	2.0 ml conical	10.2 mm	44.3 mm	1000	6823
PI0340A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	7583
PI0340A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	7518
PI0350	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	7518
PI0342A	1.5 ml amber Self-standing	10.2 mm	44 mm	1000	7583
PI0342A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	7518
PI0352	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	7518
PI0344A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	7583
PI0344A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	7518
PI0354	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	7518

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0346A	2.0 ml amber self-standing	10.2 mm	45.1 mm	1000	7583
PI0346A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	7518
PI0356	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	7518
PI0348A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	7583
PI0348A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	7518
PI0358	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	7518

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring, PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0360	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	9781
PI0361	1.5 ml Self-standing	10.2 mm	44 mm	1000	9781
PI0362	1.5 ml conical	10.2 mm	43.1 mm	1000	9781
PI0363	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	9781
PI0364	2.0 ml conical	10.2 mm	44.3 mm	1000	9781

PI0350, PI0352, PI0354, PI0356, PI0358 are available with red, blue, green, yellow, orange, purple & white color looped cap options for color coding

### Sterile Screw cap tubes with Looped Cap

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0341	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	7547
PI0343	1.5 ml Self-standing	10.2 mm	44 mm	1000	7547
PI0345	1.5 ml Conical	10.2 mm	43.1 mm	1000	7547
PI0347	2.0 ml self-standing	10.2 mm	45.1 mm	1000	7547
PI0349	2.0 ml Conical	10.2 mm	44.3 mm	1000	7547
PI0341A	0.5 ml amber Self-standing	10.2 mm	43.85 mm	1000	8394
PI0341A1	0.5 ml transparent amber Self-standing	10.2 mm	43.85 mm	1000	8394
PI0351	0.5 ml coloured Self-standing	10.2 mm	43.85 mm	1000	8394
PI0343A	1.5 ml amber Self-standing	10.2 mm	43.95 mm	1000	8394
PI0343A1	1.5 ml transparent amber Self-standing	10.2 mm	44 mm	1000	8318
PI0353	1.5 ml coloured Self-standing	10.2 mm	44 mm	1000	8394
PI0345A	1.5 ml amber conical	10.2 mm	43.1 mm	1000	8394
PI0345A1	1.5 ml transparent amber conical	10.2 mm	43.1 mm	1000	8318
PI0355	1.5 ml coloured conical	10.2 mm	43.1 mm	1000	8394

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0347A	2.0 ml amber self-standing	10.2 mm	45.1 mm	1000	8394
PI0347A1	2.0 ml transparent amber self-standing	10.2 mm	45.1 mm	1000	8318
PI0357	2.0 ml coloured self-standing	10.2 mm	45.1 mm	1000	8394
PI0349A	2.0 ml amber Conical	10.2 mm	44.3 mm	1000	8394
PI0349A1	2.0 ml transparent amber Conical	10.2 mm	44.3 mm	1000	8318
PI0359	2.0 ml coloured Conical	10.2 mm	44.3 mm	1000	8394

### Last Drop Low Retention Screw cap tubes with EPDM O-Ring,PP

Code	Capacity	Outer Diameter	Height	No./case	Price/case (₹)
PI0365	0.5 ml Self-standing	10.2 mm	43.85 mm	1000	104456
PI0366	1.5 ml Self-standing	10.2 mm	44 mm	1000	10445
PI0367	1.5 ml conical	10.2 mm	43.1 mm	1000	10445
PI0368	2.0 ml Self-standing	10.2 mm	45.1 mm	1000	10445
PI0369	2.0 ml conical	10.2 mm	44.3 mm	1000	10445

PI0351, PI0353, PI0355, PI0357 & PI0359 are available with red, blue, green, yellow, orange, purple & white color looped cap options for color coding